




MESA-50

X-ray Fluorescence Analyzer



*Speedy, Small, Simple, Smart
and Safe Elemental Analysis for you*

Touch and

What is the MESA-50 Series?


HORIBA MESA-50 Series is a portable benchtop Energy Dispersive X-ray Fluorescence (EDXRF) analyzer with user-friendly operation.

EDXRF analyzers are a fast and easy screening tool which allows elemental analysis on various sample state (solid, liquid, and powders) non-destructively.

Among them, the MESA-50 Series is the compact benchtop EDXRF analyzer with X-ray shields compliant with IEC 61010-1/JAIMA S0101-2001.

In addition to the small footprint, it is equipped with a rechargeable battery inside and it is designed as LN₂-less and pump-less design.

The MESA-50 Series enables you to perform your measurement in a variety of installation environments.



<Sample capacity>

Chamber size: 190 x 225 x 40 [W x D x H (mm)]

Maximum loading mass: up to 3 kg

Uniqueness

- Portability enabled by the internal battery
- Chamber with X-ray shields ensures your safety
- LN₂-less and pump-less optics allow for a compact design
- Recipe preset function facilitates fast and easy operation



For quick desktop analysis



For quick screening



For field work



For science education in schools



For cross validation after another analysis




For analysis in a N₂ purging glovebox*

*Please do not use it in Ar purged or He purged or other active gas purged box. It will cause discharge ignited by high voltage of X-ray generator.


Feel! MESA-50

Key Feature 1 - The Optical System



- 1 Powerful X-ray generator**
Up to 50 kV, 200 μ A, 10 W with Pd target
- 2 Switchable X-ray filters**
Enhances S/N by cutting off background
- 3 Switchable multi-collimators**
Selectable down to 1.2 mm up to 7 mm
- 4 High resolution CMOS camera**
Identifies an analysis position on a sample
- 5 Silicon drift detector**
Achieves high performance without LN₂

Key Feature 2 - The Software (V3 Software)



- 1 Auto data save**
Saves the data automatically.
It reduces the risk of losing data.
- 2 Multiple windows**
Displays multiple windows to show a XRF spectrum, periodic table, and a quantitative result at once.
- 3 Spectrum search/marker function**
Incorporates a database of the peak positions.
Users can identify the peaks in a spectrum easily even without knowledge of peak position.
- 4 Recipe function**
Allows users to save settings for measurement.
It is helpful for routine analysis.
- 5 Safety function**
Shows the status of the X-ray generator.
The interlock system inactivates a start button when the chamber door is opened.
- 6 Camera image**
Shows a measurement position on a sample,
according to the analysis spot size.



Optional Software Functions

1. RoHS Module

The RoHS module allows users to start RoHS screening and halogen screening by just one click of a button. The module offers built-in workflow which enables automatic sequential measurement without special skill and knowledge of users.

Key Functionalities

- Built-in sequences for 6 types of sample matrix classification
 - Plastic, Al alloy, Fe alloy, Cu alloy, Sn alloy, and others
- Internal parameters optimized for the 6 types of sample matrix
 - Built-in calibration curves
 - Auto thickness correction for plastic matrix
- Pass or Fail Judgement
- Auto quit mode to save time for high concentration samples
- Report template specialized for RoHS-regulated elements



2. As/Sb Module

The As/Sb module is an additional module to the RoHS module. It can determine As and Sb concentration using the built-in parameters, and export a report which includes the result of As and Sb.

Standard Accessories



Transmission X-ray film (20 sheets/pc)
Size: φ 64 mm
Material: PP
 This is a film placed at the measurement position to protect the optics from sample leakage or contamination.



Cell window (100 sheets/pc)
Size: φ 80 mm
Material: PET
 This film covers a sample cell to analyze powdery samples or liquid samples. It is also compatible with other types of cell (See the next page).



Sample cell (2 pcs)
Size: 23 x 22 mm [φ x H]
Material: PP
 This is a manually-assembled cell which can be used with a cell window to measure small fragments or liquid samples.

Check out how to assemble the cell in [the video!](#) →




3. Multilayer FPM Module

The multilayer FPM module is an optional software for thickness calculation. The software calculates the layer thickness of a multilayer sample using user-defined layer modeling. It allows users to get layer thickness results even without standard samples. The total layer thickness covering from several to tens of μm can be detected and calculated*.

*Depends on the target elements, the sample matrix, and the measurement condition.

Key Functionalities

- Two types of output: Thickness or mass-thick can be selected.
- The modeling can be defined up to 4 layers including bulk and up to 4 compositions per layer.
- It is possible to enhance the accuracy by correction using 1 or 2 standard samples.
- The layer composition can be calculated if the elements are known in each layer.



Optional Accessories



Disposable cell (48 pcs/box)

This is a pack of disposable sample cups to place a sample in liquid, cream, or powder form for analysis.

Check out how to assemble the cell in [the video!](#) ➔



Cable fixing jig (1 pc)

This is a jig to fix a cable sample on the measurement position in the sample chamber properly.





Offline license (1 pc)

This dongle of an additional offline license allows data processing on your laptop apart from the MESA-50 main unit.

Applications

1 Failure analysis / Contamination analysis



The MESA-50 Series offers you fast identification of failure or contamination happening in products, non-destructively. Thanks to the deep penetration of X-rays, you can analyze a contamination even embedded in plastic packaging as shown below. In addition, thanks to the optical camera and small spot size, the MESA-50 Series is capable of analyzing even a small defect or contamination easily.



2 Food: Salt determination in sauce



For salt determination in food, the Mohr method is typically used to determine chlorides in water. The method is based on titration and requires strong acids such as silver nitride and a skilled operator. The MESA-50 Series provides a fast and user-friendly analysis of chlorine without any strong acids and special treatment. As shown below, the results show good linearity of calibration curve plots with a regression ($R^2 > 0.99$).




3 Environment: Harmful elements in soil



The MESA-50 Series allows non-destructive analysis of soil samples without the need to destroy or dissolve the sample. Therefore, it is suitable for screening and on-site analysis. It can also analyze multiple elements, allowing the simultaneous determination of different contaminants and elements in the soil.


● Procedure



● Operation



● Analysis of harmful elements in soil



Calibration curve example

4 Steel: Coating analysis (Multilayer FPM module)



EDXRF is widely used as a non-destructive coating thickness analysis tool thanks to the analysis depth. The MESA-50 Series offers an optional software module called Multilayer FPM, which allows thickness determination of each layer of a multi-layered sample (up to 4 layers including substrate) even without the standard samples. The below example shows the analysis result of Zn-Ni alloy coating on an steel plate.



Sample: ZnNi coating (9 μm) on an steel plate

Sample	Thickness (μm)	Composition (mass%)	
		Zn	Ni
1	9.08	85.60	14.40
2	8.99	85.55	14.55
3	8.92	85.61	14.39
Average	9.00	85.59	14.41
SD	0.08	0.03	0.03

5 Electronics: RoHS screening analysis (RoHS module)



The MESA-50 Series is widely used for rapid screening to determine the presence of hazardous elements regulated by the RoHS Directive. IEC 62321-3-1 and ASTM F2617 introduce EDXRFs as a tool for determination of Cr, Br, Cd, Hg and Pb in electrical and electronic equipment, and the MESA-50 Series with the RoHS module (P.4) provides the analysis more user-friendly.




Result

Quantified: 9/20/15 3:58:58 PM



Focused Elements	Concentration	Three Sigma	Judgment	Analysis method
Pb	1194.7 ppm	57.5 ppm	NG	FFM
Cd	948.9 ppm	126.8 ppm	NG	FFM
Cr	158.3 ppm	267.6 ppm	OK	FFM
Hg	4.1 ppm	47.2 ppm	OK	FFM
Br	3.8 ppm	14.5 ppm	OK	FFM


Other Elements	Concentration	Three Sigma	Judgment	Analysis method
Sn	94.6332 wt%	0.0578 wt%	---	FFM
Ag	3.5274 wt%	0.0250 wt%	---	FFM
Cu	0.9125 wt%	0.0150 wt%	---	FFM
Mn	0.3355 wt%	0.0104 wt%	---	FFM
I	0.2794 wt%	0.0390 wt%	---	FFM
Fe	0.0183 wt%	0.0111 wt%	---	FFM
Ni	0.0027 wt%	0.0053 wt%	---	FFM
As	0.0008 wt%	0.0016 wt%	---	FFM

6 Energy: Pt loading mass determination



The MESA-50 Series allows a fast and non-destructive method for Pt catalyst loading mass determination in catalyst sheet of proton exchange membrane fuel cell. As shown on the right, several in-house samples (0.052 - 0.39 mg/cm^2) with different Pt loadings were analyzed to produce a calibration curve, which was used to determine the catalyst loading of the samples. The results showed that the curve had a good linearity and that the calculated result was consistent with the labeled value of the sample with good repeatability.





POWER

CHARGE

ERROR



Simple Operation

1 Set a sample

Sample set position



Set a sample on the analyzer.

Check the optical image




Check out the actual operation video!



Click the link!!

2 Measure and check the result



Quantitative analysis result

Measurement recipe can be saved in advance and the operator only needs to select the recipe name.

Measurement Start/Stop buttons

Outputs the reports in Excel®, Word® or PDF format

The screenshot shows a report output window with a table of results. The table has columns for 'Element', 'Concentration', and 'Thresh. Signal'. The data is the same as in the main software window.

Element	Concentration	Thresh. Signal
Zn	66.9917 wt%	0.0437 wt%
Fe	16.8550 wt%	0.0216 wt%
Ni	16.9223 wt%	0.0202 wt%

Specifications

Measurement principle	Energy dispersive X-ray fluorescence (EDXRF)
Measurable elements	Al (13) to U (92)
Sample form	Solid, liquid, powder
Chamber size	190 x 225 x 40 mm [W x D x H]
Sample observation	Color CMOS camera image
Chamber environment	Ambient
X-ray tube	Max. 50 kV, 0.2 mA
X-ray spot size	1.2 mm, 3 mm, 7 mm (automatic switching)
Filter	Four types (automatic switching)
X-ray detector	SDD (Silicon Drift Detector)
Data output	Outputs reports in Microsoft Excel®, Word® or PDF
Operating System	Microsoft Windows®11
Languages	Japanese, English, Chinese
Power supply	AC 100 to 240 ±10%, 50/60 Hz
Power consumption	0.1 kVA or less (excluding PC)
Battery	Up to 6 hours depending on operation
Weight	Approx. 12 kg
Dimensions	208 mm x 294 mm x 205 mm [W x D x H]
Options	RoHS module, As/Sb module, Multilayer FPM module, disposable cells, cable fixing jig, offline license

*MESA-50 is available in all the regions except the EU and the UK.

- The specifications, appearance or other aspects of products in this catalog are subject to change without notice.
- Please contact us with enquiries concerning further details on the products in this catalog.
- The color of the actual products may differ from the color pictured in this catalog due to printing limitations.
- It is strictly forbidden to copy the content of this catalog in part or in full.
- The screen displays shown on products in this catalog have been inserted into the photographs through compositing.
- All brand names, product names and service names in this catalog are trademarks or registered trademarks of their respective companies.
- Windows, Word and Excel are trademarks or registered trademarks of the Microsoft Corporation of the United States of America and other countries.

Other Product Line-up



MESA-50K: For larger samples



Larger Chamber



Sample size: Max. 460 × 360 × 150 mm [W × D × H]



Check out our Website!



HORIBA

<https://www.horiba.com/int/scientific/>

Japan: +81 (75) 313-8121
China: +86 (21) -6289-6060
Korea: +82 (31) 296-7911

Singapore: +65- (6) -745-8300
India: +91- (11) -4646-5000
USA: +1 (732) 494-8660

France: +33- (1) 69-74-72-00
Brazil: +55 (0) 11-2923-5400

Bulletin:HRE-3935F

Explore the future

HORIBA